
**Spencer SX**  
Barella cucchiaio

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
**Spencer SX**  
Pick Up stretcher

**CE** This appliance conforms with the Directive 93/42/CEE "Medical Devices"

Guarantee of Quality system for the production and the final control of the products certified by the notifying body TÜV SÜD Product Service GmbH

**INDEX**

General information	page 2	Operating instructions	page 6
Warnings	page 2	Maintenance and cleaning	page 9
Product description	page 5	Accessories and spare parts	page 11

First edition: 24/04/02

Rev. 6: 05/03/15

## 1. GENERAL INFORMATION

### 1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

### 1.2 Conservation of the instruction and maintenance manual

The instruction and maintenance manual must be kept together with the product, for the whole life of the device, inside a dedicated container and above all, away from any substances or liquids which could compromise perfect legibility.

### 1.3 Symbols used

Symbol	Meaning
--------	---------

General or specific warnings

See instructions for use

Lot number

Product code

The product is compliant with the specifications of the Directive 93/42/CEE

### 1.4 Servicing request

For any information regarding the correct interpretation of the instruction manual, the use, maintenance, installation and restore of the product, please contact the Spencer Customer Care Service tel. 0039 0521 541111, fax 0039 0521 541222, e-mail [service@spencer.it](mailto:service@spencer.it) or write to Spencer Italia S.r.l. – Via Provinciale, 12 - 43038 Sala Baganza (Parma) - ITALY. In order to facilitate the assistance service, please always indicate the lot number (LOT) shown on the label applied on the box or on the device.

### 1.5 Demolition

When the devices are no more suitable for being used, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste, otherwise follow the current regulations about demolition.

The disposal of the accessories and the medical device shall be performed according to the current regulations in each country.

### 1.6 Labelling

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the manufacturer, the product, the CE mark, the lot number (LOT). It must never be removed or covered.

## 2. WARNINGS

### 2.1 General warnings

- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register which will certify the eligibility of the operators to use the Spencer device has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
- Spencer Italia S.r.l. is always at your disposal to plan trainings on products.
- Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- If the instructions belong to another device and not to the device received, inform the manufacturer immediately and avoid use of the device.
- In case of any doubts about the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.

- Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
- Regularly check the appliance, carry out the prescribed maintenance and respect the average life span, as indicated by the manufacturer in this user's manual.
- Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- Use of the device in anyway other than described in this manual is forbidden.
- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- Handle with care.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments.
- When the device is being used, the assistance of qualified staff must be guaranteed.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store in a cool, dry, dark place and do not expose to direct sun.
- Store and transport device in its original packaging.
- The device not be exposed to or come into contact with any source of combustion or inflammable agents.
- Position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Attention: laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
- With reference to the D. Lgs. 24th February 1997, n. 46 emended by D. Lgs. 25/01/2010, n. 37 – Acknowledgement of Directive 93/42/CEE and 2007/47/CE, we remind both public and private operators that they are obliged to report any accident that involves any medical device to the Ministry of Health and to the Manufacture as specified and within time given by the European regulations.
- In addition, both public and private operators are obliged to inform the manufacturer of any measures that should be adopted to make the steps necessary to guarantee the safety and the health of the patients and the users o any medical device.
- As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly required to have a basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the goods final destination Country (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
- Promptly notify Spencer Italia S.r.l. regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user's manual.
- Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- The distributor or final user is aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore Spencer Italia S.r.l. expressly disclaims any responsibility and/or liability for your non-compliance with the present regulatory provisions.

## 2.2 Specific warnings

- Establish a maintenance program and periodic testing, identifying a reference employee. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in the user's manual.
- All maintenance and periodic check activities must be registered and collected together with their intervention reports (see Maintenance Register) these documents have to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
- Use only components/spare parts and/or accessories that are original or approved by Spencer Italia S.r.l. in order to carry out any operation without causing any alteration or modification to the device, otherwise we assume no responsibility for the proper functioning or damage resulting from device to the patient or the operator and warranty and will be considered void according to the compliance to the Medical Device Directive 93/42/CEE.
- Always respect the maximum load capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the load of the total weight on the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
- Never leave the patient unassisted on the device, because he may be injured.
- The device and all its components, after washing, should be allowed to dry completely before storing.
- Lubrication must be carried out after cleaning and complete drying.
- Follow the procedures approved by the Emergency Medical Service for the immobilization and transport of patients.
- Follow the procedures approved by the Emergency Medical Service for the positioning and transport of patients.
- Avoid contact with sharp objects.
- Make sure, before lifting, that the operators have a firm grip on the device.
- Avoid pulling the device on rough surfaces.
- Do not lift by crane or other mechanical lifters.
- Have practice with an empty stretcher, in order to make sure you become familiar with the manoeuvres.
- The use of the device requires at least two operators in suitable physical conditions, they must therefore be endowed with strength, balance, coordination, and common sense and they must be trained on the proper operation of the device Spencer stretcher.
- For loading techniques of the patient, for particularly heavy patients, for working on uneven ground or in special and unusual circumstances, the presence of more operators is recommended (not only 2 as expected under standard conditions).
- The maximum weight supported by each rescuer must comply with the requirements prescribed by the law of the land, in the field of Health and Safety at Work.
- Before each use, check the integrity of the belts and their hooks, as specified in the user's manual. In case of malfunction or damage that may compromise the functioning and safety of the device, patient or operator, it is necessary to replace the belts.
- Make sure the belts are properly fastened to the frame/patient board of the stretcher.
- Always immobilize the patient, using the straps supplied by the Manufacturer; lack of immobilization may cause serious damage.
- Do not operate in case the weight has not been distributed correctly.
- Use the stretcher only as described in this user's manual.
- Do not alter or modify the stretcher arbitrarily to make it fit into the ambulance: the modification may cause unforeseeable functioning and damage to the patient and operators. In any case the warranty will be void.
- Pay a lot of attention to possible obstacles (water, ice, debris, etc.) on the route of the stretcher/chair, because they could cause loss of balance of the operator and compromise the proper functioning of the device. If the path free cannot be made free of obstacles, choose an alternative path.

## 2.5 Contraindications and side effects

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

## 2.6 Physical requirements of the operators

Spencer SX is a device destined to professional use only. The rescue operators must have the following minimum requirements:

- physical capacity for operating the device
- be able to seize the device firmly with both hands
- have strong back, arms and legs for lifting, pushing and pulling the stretcher
- have a good muscular coordination

The operators must be trained in efficient, effective and safe patient transport.

Patient loading procedures for extremely heavy patients, operations in rough terrain and in particular situations more operators may be needed (not only two as in normal conditions).

**The capacities of the various operators must be considered before determining their role in the employment of the stretcher.**

**3. DESCRIPTION OF PRODUCT**

**3.1 Intended use**

The Pick Up stretcher is a device which enables the lifting and loading of a patient onto immobilising or transport devices such as spinal boards, vacuum mattresses or transport stretchers. It is essential for the loading trauma patients by two rescuers: with use of this device patient can be picked up without changing position and the trunk and limbs can maintain the supine position while loading. The structure and mechanism of the Pick Up stretcher offer maximum reduction of the movement required. The Pick Up stretcher does not offer any type of isolation. No action on behalf of the patient is expected.

**3.3 Main components**

n°	Description of component	Materials
1	Locks	Aluminium
2	Support for the lower limbs	Aluminium
3	Side blocks	Aluminium
4	Support for the body	Aluminium
5	Extensions	Aluminium


Fig. A

**3.3 Models**

**These models could be modified, with reference to codes and/or descriptions without any previous notification.**

- ST05001B Spencer SX-R 3 baldes
- ST05002B Spencer SX-R 2 baldes
- ST05040A Spencer SX-F Fixed
- ST05006A Spencer SX grey
- ST05042A Spencer SX orange
- ST05043A Spencer SX red
- ST05041A Spencer SX green
- ST05012A Spencer SX blue
- ZS00412A Spencer SX black

### 3.7 Technical data

Characteristics	Spencer SX	Spencer SX-R	Spencer SX-F
Height (mm)	70	70	70
Width (mm)	430	440	430
Minimum length (mm)	1660	1660	
1 Intermediate length (mm)	1780	1780	
1 Intermediate length (mm)	1900	1900	
Maximum length (mm)	2020	2050	1910
Folded length (mm)	1190	1190	
Weight (kg)	9,5	8	10
Load capacity (kg)	170	170	170

### 3.8 Reference standards

Reference	Title of document
MDD 93/42/CEE	European Directive about Medical Devices
MDD 2007/47/CEE	Modifications to 90/385/CEE Directive about active implants, Directive 93/42/CEE about medical devices and Directive 98/8/CE about the introduction of biocides onto the market
Legislative Decree 24/02/1997, n. 46	Application of the 93/42/CEE Directive about Medical Devices
Legislative Decree 25/01/2010, n. 35	Modifications and additions to the 20/02/97 Decree n. 46
UNI EN ISO 14971	Application of risks managing to medical devices
UNI CEI EN ISO 15223-1	Medical devices - Symbols for use in the medical device labels, labelling and information to be provided. Part 1: general requirements
UNI CEI EN 1041	Information supplied by the medical devices manufacturer
CEI EN 62366	Medical Devices - Application of the utilisation characteristics of engineering to medical devices
MEDDEV 2.4/1a-b	Guideline for the classification of medical devices
NB-MED 2.5.1/Rec 5	Technical Documentation
MEDDEV 2.7.1	Clinical Data
MEDDEV 2.12/1	Medical Devices vigilance system
UNI EN 14155	Clinical evaluation of the medical devices for human beings - Part 2: Clinical evaluation plans

### 3.9 Environmental conditions

Functioning temperature: from -10 to +40 °C

Storage temperature: from -20 to +60 °C

Relative humidity: from 5 to 95%

## 4. OPERATING INSTRUCTIONS

### 4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

Keep the original packaging for use in case of any further transport and for storage. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

### 4.2 Preparation

On receipt of the product:

- Remove the packaging and display the material so that all components are visible.
- Check that all the components/pieces on the accompanying list are present.

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage. In particular, check:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Absence of cuts, holes, tears on the structure, including the straps

- Correct fixation of all nuts, bolts and screws
- Correct fixation of straps
- Correct fastening of straps
- State of use (moving parts, belts)
- Integrity of components
- Lubrication of moving parts
- There are seat belts for the immobilization of the patient and they are intact and functioning
- No piping or metal sheet present bends or cracks
- Test the hooks closing them and opening them in succession
- Check cleanliness of the hooks and the absence of extraneous parts to the mechanisms

If the above conditions are met, the device may be considered ready for use; otherwise you must immediately remove the device from service and contact the Manufacturer.

### 4.3 Functioning

#### 4.3.1 Loading the patient

**A correct medical evaluation should be made before the patient is moved, lifted or transported.**

**Before loading the patient the stretcher must be positioned as close as possible to the patient.**

**Before proceeding with transport, it is essential that the patient is fastened to the Pick Up stretcher using the belts supplied which will offer stability and safety.**

**The Pick Up stretcher is not indicated for immobilising or transporting for long periods / distance**

- Start with the initial immobilising of the cervical area as prescribed by your standard area protocols.
- Position the Pick Up stretcher in the closed position as close as possible to the patient. The lengthening of the top support profile should be reduced to a minimum and it should correspond to the head of the patient (fig. B).
- To lengthen the profile - turn the side blocks.
- Lengthen the supporting profile for the lower limbs so that it corresponds with the patient heels and turn the side locks until you hear the locks close (fig. C).
- After checking that the locks are correctly closed and the length is correctly adjusted open the SPENCER SX RED KEY vertically and position them on the sides of the patient (fig D).
- The operators should position themselves on the sides of the patient in the trunk area. Move each half of the device and position them under the patient. A minimum movement of a few centimetres in the shoulder and the hip area can be made in order to assist this operation (fig E).
- In order to lock - proceed first with the head lock and then the area of the feet making sure that the red lock is pushed in so that the two halves lock together correctly (fig. F).

**Before moving the device with the patient loaded, check that the two halves of the stretcher are correctly attached to each other both at the foot end and at the head end. Check also that the length adjusting mechanisms are correctly blocked.**

**The correct alignment of the cervical rachis and the vertebral column must be maintained.**

**The insertion of the two halves under the patient could be impeded or hampered by the ground underneath the patient or by clothes.**

**During the closure of the two halves of the device, check that no parts of the patient or of his clothes are caught up.**


Fig. B


Fig. C


Fig. D


Fig. E


Fig. F

Fig. G

**During the closure of the two halves of the device, check that no parts of the patient or of his clothes are caught up.**

- Cover the patient with an isothermal sheet which should also be tucked around the patient's body.
- The patient should be fixed to the device bearing in mind that the support is only for temporary use previous to transfer to the correct type of support (spinal board or vacuum mattress).

**Lifting and transport of the device with the patient loaded must be carried out by at least two rescuers bearing in mind the maximum load that they are permitted to lift and the physical structure and condition of the operator. If the Pick Up stretcher is lifted from the head / foot area the resistance to flexion of the stretcher decreases; for this reason we advise to lift the Pick Up stretcher from the sides (fig. H).**

**Owing to the structural characteristics of this kind of stretcher we do not advise to keep a trauma patient on the Pick Up stretcher for any more time than essential or during ambulance transport and we advice to move and transport the patient on a spinal board or a vacuum mattress.**


Fig. H

#### 4.3.2 End of use

Once the transport of the patient has terminated, carry out the procedures described in paragraph regarding maintenance and cleaning before putting the device into storage.

The Pick Up stretcher should be reduced to its minimal size for storage by folding the upper half on top of the lower half.

- Open both of the side locks and pull out the two halves until the joints are exposed (fig. I) and the Pick Up stretcher reaches its maximum extension before folding the stretcher (fig. L).

If the joints are not completely extracted the stretchers could be irreparably damaged.

- Place the Pick Up stretcher into its transport bag.

If the Pick Up stretcher is stored with its forks (vertical and horizontal) it must be shortened to its minimal length before being fixed to the fork.


Fig. I


Fig. L

#### 4.4 Troubleshooting

PROBLEM	CAUSE	REMEDY
One or all the locks do not open/close (without patient loaded)	The spring in the block is damaged	Put the device out of service and contact the Assistance Center
	Blocks are dirty	Accurately clean blocks
The extension for the lower limbs does not lengthen/shorten after unlocking the lock	Blocks are dirty	Accurately clean blocks
	Deformed profiles	Put the device out of service and contact the Assistance Center
The extension for the lower limbs does not bend up to close for storage	The joint has not been completely exposed	Lengthen the stretcher more until the joints are completely exposed.
The extension for the lower limbs does not block at the desired length	There are no holes to block the hinge	Lengthen the stretcher more until it blocks
	Blocking system is broken	Put the device out of service and contact the Assistance Center
The supporting profiles turn around the borders.	Fixing rivets are broken	Put the device out of service and contact the Assistance Center
Damage to the structure	Improper use	Put the device out of service and contact the Assistance Center

### 5. MAINTENANCE AND CLEANING

#### 5.1 Cleaning

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.

**The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures.**

The exposed metal parts are usually treated and/or painted in order to increase their resistance.

The board has been made out of incontaminable material, in order to increase hygiene and easy cleaning.

Clean the exposed parts with water and delicate soap, then dry with a soft cloth. In order to obtain a shine effect, it is possible to use car waxes and creams.

Do not clean with high pressure water; this will damage the joints and the lubricated parts.

If the stretcher is not cleaned regularly, this may cause risks in terms of cross-contamination.

We recommend the use of the polishing detergent Spencer STX 99.

Rinse thoroughly with warm water making sure that you have removed all traces of detergent, which could degrade or compromise the integrity and durability of the device. The use of high pressure water should be avoided. Water penetrates the joints and removes the oil, creating the risk of corrosion of components.

Allow to dry thoroughly before storing. Drying after washing or after use in wet environments must be natural and not forced, do not use flames or other sources of direct heat.

**In the presence of blood, oxidize it before to washing the device with water.**

## 5.2 Maintenance

### 5.2.1 Precautionary maintenance

The person who carries out the precautionary maintenance of the appliance (user in person, Manufacturer/supplier or a third party) has to guarantee the following basic requirements:

- Technical knowledge of the appliance and of the periodic maintenance procedures as described in these instructions.
- Specific qualifications and training in the maintenance operations of the appliance in question.
- The use of components/replacement parts/accessories that are either original or approved by the supplier, in such a way that each operation causes no alteration or modification to the appliance.
- Possession of the checklist of operations carried out on the appliance.
- Guarantee complete adherence to the instructions of the Directive 93/42/CEE which includes also the obligation towards the Manufacturer to maintain post sales records and traceability of the appliance if requested.

**The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures.**

Checks to be carried out before and after each use, and at least every month, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Absence of cuts, holes, tears on the structure, including the straps
- Correct fixation of all nuts, bolts and screws
- Correct fixation of straps
- Correct fastening of straps
- State of use (moving parts, belts)
- Integrity of components
- Lubrication of moving parts
- There are seat belts for the immobilization of the patient and they are intact and functioning
- No piping or metal sheet present bends or cracks
- Test the hooks closing them and opening them in succession
- Check cleanliness of the hooks and the absence of extraneous parts to the mechanisms

**The inspection frequency is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.** Please note that you must do the cleaning as described in paragraph 5.1 and verify functionality before and after each use. Spencer Italia S.r.l. declines any responsibility for the proper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance warranty and will void the compliance to the Medical Device Directive 93/42/CEE.

The person responsible for routine maintenance must identify damaged/worn parts, but the replacement or restoration of them can only be carried out by the manufacturer or by an authorized service centre.

For other replacement/repair activities contact the Manufacturer or an authorized centre.

Use only accessories/original spare parts approved by Spencer Italia S.r.l., otherwise we will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres. Warranty will be considered void in compliance with the Medical Device Directive 93/42/EEC.

### 5.2.2 Periodic maintenance

Planned interventions by the Manufacturer or authorized center are not required, but it is prescribed to make cleaning and checking indicated in the specific sections "Cleaning" and "Precautionary Maintenance".

### 5.2.3 Special servicing

**Only the Manufacturer or centres with written authorisation are authorised to complete any special servicing operations.**

For any operations that are not carried out directly by the Manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device, if used as indicated in the following instruction manual, has an average life span of 5 years.

Spencer Italia S.r.l. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the Manufacturer or by one of the Manufacturer's Authorised Service centres, making void both the guarantee and the conformity to the Medical Devices Directive 93/42/CEE.

## **6. ACCESSORIES AND SPARE PARTS**

### **6.3 Accessories**

ST05009A	SX Bag - Carrying case for Pick Up stretchers, orange
ST05010A	SX Bag - Carrying case for Pick Up stretchers, blue
ST05016A	SXO - Horizontal fixing system in aluminium for Pick Up stretcher
ST05017A	SXV - Vertical fixing system in aluminium for Pick Up stretcher
ST05018B	SXV – Vertical fixing system in aluminium for Pick Up stretcher 20G
ST05019B	SXO – Horizontal fixing system in aluminium for Pick Up stretcher 20G
ST05050B	Ranger SX – Foldable trolley for Pick Up stretchers' transport
ST00598A	STX 598 – Two pieces belt with plastic hook and orange belt
ST00592A	STX 592 - Two pieces belt with metal hook and orange belt
ST00500C	Strap Up – Belt cover, blue

### **6.4 Spare parts**

ST05048A	SX-R Red Key/Z5000 - Plastic joint for replacement
----------	----------------------------------------------------

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© Copyright Spencer Italia S.r.l.